

GE Healthcare
Life Sciences

ÄKTAexplorer™, ÄKTApurifier™ and ÄKTAmicro™

Operating Instructions

Original instructions



Table of Contents

| | | |
|----------|--|-----------|
| 1 | Introduction | 5 |
| 1.1 | Important user information | 6 |
| 1.2 | Regulatory information | 7 |
| 1.3 | Instrument | 10 |
| 1.4 | Control software | 14 |
| 1.5 | User documentation | 14 |
| 2 | Safety instructions | 16 |
| 2.1 | Safety precautions | 16 |
| 2.2 | Labels | 24 |
| 2.3 | Emergency procedures | 26 |
| 2.4 | Recycling information | 27 |
| 2.5 | Declaration of Hazardous Substances (DoHS) | 29 |
| 3 | Installation | 32 |
| 3.1 | Site requirements | 32 |
| 3.2 | Transport | 33 |
| 3.3 | Unpacking | 34 |
| 3.4 | Assembly | 34 |
| 3.5 | Connections | 35 |
| 3.6 | Spare parts and accessories | 37 |
| 4 | Operation | 38 |
| 4.1 | Operation overview | 38 |
| 4.2 | Starting the instrument | 38 |
| 4.3 | Starting the control system | 39 |
| 4.4 | Preparations before start | 40 |
| 4.5 | Setting up a run | 43 |
| 4.6 | Performing a run | 47 |
| 4.7 | Procedures after a run | 50 |
| 5 | Maintenance | 51 |
| 5.1 | General | 51 |
| 5.2 | User maintenance schedule | 52 |
| 5.3 | Cleaning | 54 |
| 5.4 | Component maintenance | 56 |
| 5.5 | Disassembly and assembly of components and consumables | 56 |
| 5.6 | Replacement of fuses | 56 |
| 5.7 | Calibration | 57 |
| 5.8 | Storage | 58 |
| 6 | Troubleshooting | 59 |
| 6.1 | UV curve problems | 59 |
| 6.2 | Conductivity curve problems | 60 |

Table of Contents

| | | |
|----------|--|-----------|
| 6.3 | pH curve problems | 62 |
| 6.4 | Pressure curve problems | 65 |
| 7 | Reference information | 66 |
| 7.1 | Specifications | 66 |
| 7.2 | Chemical resistance | 66 |
| 7.3 | System recommendations | 70 |
| 7.4 | Health and Safety Declaration Form | 71 |
| 7.5 | Ordering information | 73 |
| A | Connection diagram - Liquid flow path | 74 |
| B | Tubing | 78 |

1 Introduction

Purpose of the Operating Instructions

The *Operating Instructions* provide you with the instructions needed to handle ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro in a safe way.

Prerequisites

In order to operate ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro as is intended, the following pre-requisites must be fulfilled:

- The user should have a general understanding of how a PC and the Microsoft™ Windows™ operating system works.
 - The user must understand the concepts of liquid chromatography.
 - The user must read and understand the Safety Instructions in this manual.
 - ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro and software should be installed, configured and calibrated according to these Operating Instructions.
-

About this chapter

This chapter contains important user information, a description of the intended use of ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro, regulatory information, list of associated documentation, definitions of safety notices and so on.

1.1 Important user information

Read this before operating ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro



All users must read the entire *Operating Instructions* before installing, operating or maintaining ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro.

Always keep the *Operating Instructions* at hand when operating ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro.

Do not operate ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro in any other way than described in the user documentation. If you do, you may be exposed to hazards that can lead to personal injury and you may cause damage to the equipment.

Intended use

ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro are liquid chromatography systems intended for protein purification within method development and drug discovery. The systems can be used to screen for optimal choice of columns, media and running parameters to purify selected proteins.

ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro are intended for research use only, and shall not be used in any clinical procedures, or for diagnostic purposes.

Safety notices

This user documentation contains WARNINGS, CAUTIONS and NOTICES concerning the safe use of the product. See definitions below.

Warnings



WARNING

WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury. It is important not to proceed until all stated conditions are met and clearly understood.

Cautions



CAUTION

CAUTION indicates a hazardous situation which, if not avoided, could result in minor or moderate injury. It is important not to proceed until all stated conditions are met and clearly understood.

Notices



NOTICE

NOTICE indicates instructions that must be followed to avoid damage to the product or other equipment.

Notes and tips

Note: *A note is used to indicate information that is important for trouble-free and optimal use of the product.*

Tip: *A tip contains useful information that can improve or optimize your procedures.*

Typographical conventions

Software items are identified in the text by ***bold italic*** text. A colon separates menu levels, thus ***File:Open*** refers to the ***Open*** command in the ***File*** menu.

Hardware items are identified in the text by **bold** text (e.g., **Power** switch).

1.2 Regulatory information

In this section

This section describes the directives and standards that are fulfilled by ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro.

Manufacturing information

The table below summarizes the required manufacturing information. For further information, see the EC Declaration of Conformity document.

| Requirement | Content |
|----------------------------------|--|
| Name and address of manufacturer | GE Healthcare Bio-Sciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden |

CE conformity

This product complies with the European directives listed in the table, by fulfilling the corresponding harmonized standards.

A copy of the EC Declaration of Conformity is available on request.

| Directive | Title |
|-------------|---|
| 2006/42/EC | Machinery Directive (MD) |
| 2006/95/EC | Low Voltage Directive (LVD) |
| 2004/108/EC | Electromagnetic Compatibility (EMC) Directive |

CE marking



The CE marking and the corresponding Declaration of Conformity is valid for the instrument when it is:

- used as a stand-alone unit, or
 - connected to other CE marked instruments, or
 - connected to other products recommended or described in the user documentation, and
 - used in the same state as it was delivered from GE Healthcare, except for alterations described in the user documentation.
-

International standards

This product fulfills the requirements of the following standards:

| Standard | Description | Notes |
|---|---|--|
| EN/IEC 61010-1, UL 61010-1, CAN/CSA-C22.2 No. 61010-1 | Safety requirements for electrical equipment for measurement, control, and laboratory use | EN standard is harmonized with EU directive 2006/95/EC |
| EN 61326-1 | Electrical equipment for measurement, control and laboratory use - EMC requirements | EN standard is harmonized with EU directive 2004/108/EC |
| EN ISO 12100 | Safety of machinery. General principles for design. Risk assessment and risk reduction. | EN ISO standard is harmonized with EU directive 2006/42/EC |

Regulatory compliance of connected equipment

Any equipment connected to ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro should meet the safety requirements of EN 61010-1/IEC 61010-1, or relevant harmonized standards. Within the EU, connected equipment must be CE marked.

Environmental conformity

| Requirement | Title |
|------------------------------|--|
| 2011/65/EU | Restriction of Hazardous Substances (RoHS) Directive |
| 2012/19/EU | Waste Electrical and Electronic Equipment (WEEE) Directive |
| Regulation (EC) No 1907/2006 | Registration, Evaluation, Authorization and restriction of Chemicals (REACH) |
| ACPEIP | Administration on the Control of Pollution Caused by Electronic Information Products, China Restriction of Hazardous Substances (RoHS) |

1.3 Instrument

Product description

ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro are high pressure liquid chromatography systems for use in laboratory scale production of biomolecules.





Figure 1.1: The main parts of the instruments. ÄKTAexplorer is shown as an example.

| Part | Function |
|------|-------------------------------|
| 1 | Monitor pH/C-900 |
| 2 | Monitor UV-900 |
| 3 | Pump P-900 (P-901 alt. P-903) |
| 4 | Power switch |
| 5 | Switch valve (SV-903) |
| 6 | Conductivity cell |
| 7 | UV cell |

| Part | Function |
|------|---|
| 8 | Mixer M-925 |
| 9 | Sample pump P-960 (This is only standard in some versions) |
| 10 | Valve door |
| 11 | Column selection valve, V3 (PV-908) |
| 12 | Sample valve, V5 (PV-908) |
| 13 | Injection valve, V1 (INV-907) |
| 14 | Column selection valve, V2 (PV-908) |
| 15 | Flow direction valve, V7 (INV-907). This is only standard for ÄKTAEplorer 100 and ÄKTAEplorer 100 Air |
| 16 | Outlet valve, V4 (PV-908) |
| 17 | Buffer valve, V6 (INV-908) |
| 18 | Box 900 |

Detailed information on the components included in each system can be found in their respective User Manual.

Electrical and communication connections



| No. | Description | No. | Description |
|-----|-------------------------------|-----|-------------------------|
| 1 | Fraction collector (optional) | 5 | CU-950 |
| 2 | Mains supply socket | 6 | UniNet-1-cable |
| 3 | UniNet-1-cable | 7 | USB cable (to computer) |
| 4 | Power cord | 8 | Power converter |

Basic flow path




Figure 1.2: Basic flow path. ÄKTAexplorer is shown as an example.

| Step | Part | Description |
|------|---------------|--|
| 1 | B | Buffers from container pass through Buffer valve V6 and/or the switch valve. |
| 2 | S | Sample pass through sample valve V5 which selects a sample depending on the setting in the control software. |
| 3 | pA | Pump A pumps buffer through the system. |
| 4 | M | The buffers pass through Mixer M where they are mixed. |
| 5 | V1, p1 | This is where the sample is added to the flow path. The sample can be added manually with a suitable syringe or pumped in via the sample pump p1 from a sample chosen via sample valve V5 or an autosampler. |
| 6 | V7 | The flow direction valve V7 is optional in some systems and used to select the direction of the flow through column. |
| 7 | V2, V3 | Column selection valves V2 and V3 direct flow through a specified column. |
| 8 | V7, UV, C, pH | Liquid returns to flow direction valve V7 and redirects flow to outlet valve via pH, UV and Conductivity monitors. |
| 9 | V4 | The outlet valve directs the flow either to waste, to fraction collection containers or to a fraction collector device such as Frac-950. |

The flow path is likely to include air sensors, flow restrictors, online filters, sample loops and so on which vary in number and application depending on the system and its strategy.

For more details on liquid flow path, see *Appendix A Connection diagram - Liquid flow path, on page 74*

1.4 Control software

UNICORN™ control software

UNICORN is a complete software for control and supervision of ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro. The software runs under Microsoft™ Windows operating system.

UNICORN is supplied with a method wizard which provides easy creation of methods for purification.

For more information about UNICORN control system, see the UNICORN user manuals supplied.

1.5 User documentation

In addition to these *Operating Instructions*, the documentation package supplied with ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro also includes product documentation binders containing detailed specifications and traceability documents.

The most important documents in the document package with regard to technical aspects of ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro are:

System-specific documentation

| User documentation | Content |
|--|---|
| <i>ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro Operating Instructions</i> | All instructions needed to operate the instrument in a safe way, including brief system description, installation, and maintenance. |
| <i>ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals</i> | Detailed system description. Comprehensive user instructions, method creation, operation, advanced maintenance and troubleshooting. |
| <i>ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro Installation Guides</i> | Instructions for installation and installation test. |

| User documentation | Content |
|---|---|
| EC Declaration of Conformity for ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro | Document whereby the manufacturer ensures that the product satisfies and is in conformity with the essential requirements of the applicable directives. |

Software documentation

Together with each system, the following software documentation is supplied providing additional information that applies to ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro, independent of the specific configuration:

| Document | Purpose/Contents |
|--------------------------------|---|
| <i>UNICORN™ manual package</i> | <ul style="list-style-type: none"> The manuals contain detailed instructions on how to administer UNICORN, work with methods, perform runs and evaluate results. The Online help contains dialog descriptions for UNICORN. The Online help is accessed from the Help menu. |

Component documentation

Documentation for components produced both by GE Healthcare and by a third-party are, if existent, also included in the document package.

2 Safety instructions

About this chapter

This chapter describes safety compliance, safety labels, general safety precautions, emergency procedures, power failure and recycling of ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro.

2.1 Safety precautions

Introduction

The ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instruments are powered by mains voltage and handles pressurized liquids that may be hazardous. Before installing, operating or maintaining the system, you must be aware of the hazards described in this manual. **Follow the instructions provided to avoid personal injury or damage to the equipment.**

The safety precautions in this section are grouped into the following categories:

- General precautions
 - Using flammable liquids
 - Personal protection
 - Installing and moving the instrument
 - System operation
 - Maintenance
-

General precautions

Always follow these General precautions to avoid injury when using the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instrument.



WARNING

Do not operate ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro in any other way than described in the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro and UNICORN manuals. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.



WARNING

Operation and user maintenance of the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instrument should be performed by properly trained personnel only.



WARNING

Before connecting a column to the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instrument, read the instructions for use of the column. To avoid exposing the column to excessive pressure, make sure that the pressure limit is set to the specified maximum pressure of the column.



WARNING

Do not use any accessories not supplied or recommended by GE Healthcare.



WARNING

Do not use ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro if it is not working properly, or if it has suffered any damage, for example:

- damage to the power cord or its plug
- damage caused by dropping the equipment
- damage caused by splashing liquid onto it

2 Safety instructions

2.1 Safety precautions



CAUTION

Waste tubes and containers must be secured and sealed to prevent accidental spillage.



CAUTION

Make sure that the waste container is dimensioned for maximum possible volume when the instrument is left unattended.



NOTICE

Avoid condensation by letting the unit equilibrate to ambient temperature.

Using flammable liquids

When using flammable liquids with the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instrument, follow these precautions to avoid any risk of fire or explosion.



WARNING

Fire Hazard. Before starting the system, make sure that there is no leakage.



WARNING

A fume hood or similar ventilation system shall be installed when flammable or noxious substances are used.

Personal protection



WARNING

Always use appropriate Personal Protective Equipment (PPE) during operation and maintenance of ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro system.



WARNING

When using hazardous chemical and biological agents, take all suitable protective measures, such as wearing protective glasses and gloves resistant to the substances used. Follow local and/or national regulations for safe operation and maintenance of ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro.



WARNING

Spread of biological agents. The operator has to take all necessary actions to avoid spreading hazardous biological agents in the vicinity of the instrument. The facility should comply with the national code of practice for biosafety.



WARNING

High pressure. ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro operates under high pressure. Wear protective glasses and other required Personal Protective Equipment (PPE) at all times.



WARNING

Personal Protective Equipment (PPE). Whenever packing, unpacking, transporting or moving ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro, wear:

- Protective footwear, preferably with steel lining.
- Working gloves, protecting against sharp edges.
- Protective glasses.

Installing and moving the instrument



WARNING

Supply voltage. Make sure that the supply voltage at the wall outlet corresponds to the marking on the instrument, before connecting the power cord.



WARNING

Protective ground. ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro must always be connected to a grounded power outlet.



WARNING

Power cord. Only use power cords with approved plugs delivered or approved by GE Healthcare.



WARNING

Access to power switch and power cord with plug. Do not block access to the power switch and power cord. The power switch must always be easy to access. The power cord with plug must always be easy to disconnect.



WARNING

Installing the computer. The computer should be installed and used according to the instructions provided by the manufacturer of the computer.



CAUTION

Heavy object. Use suitable lifting equipment when moving the systems. Three people are required to lift the system safely.



NOTICE

Any computer used with the equipment shall comply with IEC 60950 and be installed and used according to the manufacturer's instructions.



NOTICE

Disconnect power. To prevent equipment damage, always disconnect power from the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instrument before an instrument module is removed or installed, or a cable is connected or disconnected.



NOTICE

ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro shall be installed and prepared by GE Healthcare personnel or third party authorized by GE Healthcare.

System operation



WARNING

Hazardous chemicals during run. When using hazardous chemicals, run **System CIP** and **Column CIP** to flush the entire system tubing with distilled water, before service and maintenance.



WARNING

Hazardous biological agents during run. When using hazardous biological agents, run **System CIP** and **Column CIP** to flush the entire system tubing with bacteriostatic solution (e.g. NaOH) followed by a neutral buffer and finally distilled water, before service and maintenance.

2 Safety instructions

2.1 Safety precautions



CAUTION

Hazardous chemicals in UV flow cell. Make sure that the entire flow cell has been flushed thoroughly with bacteriostatic solution, for example NaOH, and distilled water, before service and maintenance.



NOTICE

Do not run **Column CIP** if using silica based packing material or RPC columns. Remove the column from the system during CIP.

Maintenance



WARNING

Electrical shock hazard. All repairs should be done by service personnel authorized by GE Healthcare. Do not open any covers or replace parts unless specifically stated in the user documentation.



WARNING

Disconnect power. Always disconnect power from the instrument before replacing any component on the instrument, unless stated otherwise in the user documentation.



WARNING

Hazardous chemicals during maintenance. When using hazardous chemicals for system or column cleaning, wash the system or columns with a neutral solution in the last phase or step.



WARNING

Do not perform any type of maintenance work while the system is powered electrically or when the piping system is pressurized. Note that the piping system can be pressurized even when the system is closed down.



WARNING

Only spare parts and accessories that are approved or supplied by GE Healthcare may be used for maintaining or servicing ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro.



WARNING

Make sure that the piping system is completely leakage free before performing any CIP on the system.



WARNING

NaOH is corrosive and therefore dangerous to health. When using hazardous chemicals, avoid spillage and wear protective glasses and other suitable Personal Protective Equipment (PPE).



WARNING

After assembly, the piping system must be tested for leakage at maximum pressure for continued protection against injury risks due to fluid jets, burst pipes or explosive atmosphere.



WARNING

Before disassembly, check that there is no pressure in the piping system.

2 Safety instructions

2.1 Safety precautions



WARNING

Disconnect power. Always disconnect power from the instrument before replacing fuses.



WARNING

Decontaminate the equipment before decommissioning to ensure that hazardous residues are removed.



NOTICE

Cleaning. Keep the instrument dry and clean. Wipe regularly with a soft damp tissue and, if necessary, a mild cleaning agent. Let the instrument dry completely before use.


2.2 Labels

In this section

This section describes the instrument labels and labels concerning hazardous substances that are attached to the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instruments. For information about marking of the computer equipment, refer to the manufacturer's instructions.

Labels on the instrument



The illustration below shows an example of the identification label that is attached to the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instruments.



Symbols used in instrument labels

| Label | Meaning |
|-------|---|
| | Warning! Read the user documentation before using the system. Do not open any covers or replace parts unless specifically stated in the user documentation. |
| | The system complies with the requirements for electromagnetic compliance (EMC) in Australia and New Zealand. |
| | The system complies with applicable European directives. |
| | This symbol indicates that ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro has been certified by a Nationally Recognized Testing Laboratory (NRTL). NRTL means an organization, which is recognized by the US Occupational Safety and Health Administration (OSHA) as meeting the legal requirements of Title 29 of the Code of Federal Regulations (29 CFR), Part 1910.7. |

Labels concerning hazardous substances

| Label | Meaning |
|---|---|
|  | This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of equipment. |
|  | This symbol indicates that the product contains hazardous materials in excess of the limits established by the Chinese standard SJ/T11363-2006 Requirements for Concentration Limits for Certain Hazardous Substances in Electronics. |


2.3 Emergency procedures

In this section

This section describes how to do an emergency shutdown of the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro systems. The section also describes the result in the event of power failure.

Emergency shutdown

In an emergency situation, do as follows to stop the run:

| Step | Action |
|------|--|
| 1 | To pause the run from UNICORN, click the Pause button in System Control .  |
| 2 | If required, switch off power to the instrument by pressing the Main power switch to the 0 position. |

Power failure

The result of a power failure depends on which unit that is affected.

| Power failure to... | will result in... |
|--|--|
| ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro | <ul style="list-style-type: none">• The run is interrupted immediately, in an undefined state• The data collected up to the time of the power failure is available in UNICORN |
| Computer | <ul style="list-style-type: none">• The UNICORN computer shuts down in an undefined state• The run continues, but data cannot be saved in UNICORN. |

2.4 Recycling information

Decontamination

ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro shall be decontaminated before decommissioning and all local regulations shall be followed with regard to scrapping of the equipment.

Disposal, general instructions

When taking ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro out of service, the different materials must be separated and recycled according to national and local environmental regulations.

Recycling of hazardous substances

ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instruments contain hazardous substances. Detailed information is available from your GE Healthcare representative.

Disposal of electrical components

Waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of equipment.



2.5 Declaration of Hazardous Substances (DoHS)


Introduction

The following product pollution control information is provided according to SJ/T11364-2006 Marking for Control of Pollution caused by Electronic Information Products.

根据SJ/T11364-2006《电子信息产品污染控制标识要求》特提供如下有关污染控制方面的信息


Symbols used in pollution control label

电子信息产品污染控制标志说明

| Label | Meaning |
|---|---|
|  | <p>This symbol indicates the product contains hazardous materials in excess of the limits established by the Chinese standard SJ/T11363-2006 Requirements for Concentration Limits for Certain Hazardous Substances in Electronic Information Products. The number in the symbol is the Environment-friendly Use Period (EFUP), which indicates the period during which the toxic or hazardous substances or elements contained in electronic information products will not leak or mutate under normal operating conditions so that the use of such electronic information products will not result in any severe environmental pollution, any bodily injury or damage to any assets. The unit of the period is "Year".</p> <p>In order to maintain the declared EFUP, the product shall be operated normally according to the instructions and environmental conditions as defined in the product manual, and periodic maintenance schedules specified in Product Maintenance Procedures shall be followed strictly.</p> <p>Consumables or certain parts may have their own label with an EFUP value less than the product. Periodic replacement of those consumables or parts to maintain the declared EFUP shall be done in accordance with the Product Maintenance Procedures.</p> <p>This product must not be disposed of as unsorted municipal waste, and must be collected separately and handled properly after decommissioning.</p> |

2 Safety instructions

2.5 Declaration of Hazardous Substances (DoHS)

| Label | Meaning |
|---|--|
|  | <p>该标志表明本产品含有超过SJ/T11363-2006《电子信息产品中有毒有害物质的限量要求》中限量的有毒有害物质。标志中的数字为本产品的环保使用期，表明本产品在正常使用的条件下，有毒有害物质不会发生外泄或突变，用户使用本产品不会对环境造成严重污染或对其人身、财产造成严重损害的期限。单位为年。</p> <p>为保证所申明的环保使用期限，应按产品手册中所规定的环境条件和方法进行正常使用，并严格遵守产品维修手册中规定的期维修和保养要求。</p> <p>产品中的消耗件和某些零部件可能有其单独的环保使用期限标志，并且其环保使用期限有可能比整个产品本身的环保使用期限短。应到期按产品维修程序更换那些消耗件和零部件，以保证所申明的整个产品的环保使用期限。</p> <p>本产品在使用寿命结束时不可作为普通生活垃圾处理，应被单独收集妥善处理</p> |

List of hazardous substances and their concentrations

产品中有毒有害物质或元素的名称及含量

Indication for each major part if substance exceeds limit

| Value | Meaning |
|-------|--|
| O | <p>Indicates that this toxic or hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in SJ/T11363-2006.</p> <p>表示该有毒有害物质在该部件所有均质材料中的含量均在SJ/T11363-2006 标准规定的限量要求以下</p> |
| X | <p>Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in SJ/T11363-2006.</p> <ul style="list-style-type: none">Data listed in the table represents best information available at the time of publication <p>表示该有毒有害物质至少在该部件的某一均质材料中的含量超出 SJ/T11363-2006 标准规定的限量要求</p> <ul style="list-style-type: none">此表所列数据为发布时所能获得的最佳信息 |

List of hazardous substances

| Component name 部件名称 | Hazardous substance 有毒有害物质或元素 | | | | | |
|--|----------------------------------|---------|---------|-------------|-------------|---------------|
| | Pb 铅 | Hg 汞 | Cd 镉 | Cr6+ 六价铬 | PBB 多溴联苯 | PBDE 多溴二苯醚 |
| ÄKTExplorer10S, 18-1145-05 ¹ | X | 0 | 0 | 0 | 0 | 0 |
| ÄKTExplorer100, 18-1112-41 ¹ | X | 0 | 0 | 0 | 0 | 0 |
| ÄKTExplorer100 AIR, 18-1403-00 ¹ | X | 0 | 0 | 0 | 0 | 0 |
| ÄKTExplorer100 Crystal Custom, 18-1900-71 ¹ | X | 0 | 0 | 0 | 0 | 0 |
| ÄKTApurifier10, 28-4062-64 ¹ | X | 0 | 0 | 0 | 0 | 0 |
| ÄKTApurifier100, 28-4062-66 ¹ | X | 0 | 0 | 0 | 0 | 0 |
| ÄKTApurifierUPC 10, 28-4062-68 ¹ | X | X | 0 | 0 | 0 | 0 |
| ÄKTApurifierUPC 100, 28-4062-71 ¹ | X | X | 0 | 0 | 0 | 0 |
| ÄKTAmicro, 28-9483-03 ¹ | X | X | 0 | 0 | 0 | 0 |

¹ The product has not been tested as per the Chinese standard SJ/T11363-2006 *Requirements for Concentration Limits for Certain Hazardous Substances in Electronic Information Product*.

3 Installation



NOTICE

ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro shall be installed and prepared by GE Healthcare personnel or third party authorized by GE Healthcare.



NOTICE

Any computer used with the equipment shall comply with IEC 60950 and be installed and used according to the manufacturer's instructions.

ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro are delivered in protective packing material and shall be unpacked with great care.

Any equipment connected to ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro must fulfill applicable standards and local regulations.

For detailed information on Installation, see ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals.

3.1 Site requirements

| Parameter | Requirement |
|---------------------|--|
| Electrical power | 100-240 V, 50-60 Hz |
| Ambient temperature | 4°C to 40°C |
| Placement | Stable laboratory bench min. 200 × 80 cm |
| Humidity | 20% to 95%, non-condensing |

3.2 Transport



WARNING

Personal Protective Equipment (PPE). Whenever packing, unpacking, transporting or moving ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro, wear:

- Protective footwear, preferably with steel lining.
- Working gloves, protecting against sharp edges.
- Protective glasses.

The equipment weights are specified in the table below. Each system requires at least three people to lift and move it unless a suitable lifting device is used.

| Instrument | Weight |
|--------------|--------|
| ÄKTAexplorer | 66 kg |
| ÄKTApurifier | 41 kg |
| ÄKTAmicro | 55 kg |

The instrument can be transported on a trolley or a suitable lifting device capable of supporting the weight of the instrument.



NOTICE

Lift the instrument in the upright position. Do not use the front panel bar as a lifting handle.

Before moving the system:

- disconnect all cables and tubing connected to peripheral components and liquid containers.
- remove all items from the top of the instrument.
- close the valve door completely (only for ÄKTAexplorer).
- grasp the system firmly by placing the fingers in the gap between the swivel platform and the base of the main unit and lift.

For more information on transport, see ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals.

3.3 Unpacking

Check for damage

Check the equipment for damage before starting assembly and installation. There are no loose parts in the transport box. All parts are either mounted on the system or located in the accessory kit box. If any damage is found, document the damage, and contact your local GE Healthcare representative.

Unpack the system

Remove straps and packing material. Then stand equipment upright on swivel foot before starting installation.

3.4 Assembly

The following parts must be added to the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro instruments before they can be used:

- pH electrode (optional)
- Waste tube
- CU-950 Control unit between unit and computer
- Various buffer and sample bottles

Installing the pH electrode

Install the pH electrode in the flow cell according to the picture below.




Figure 3.1: Main parts of the pH electrode and holder.

| No. | Description | No. | Description |
|-----|-------------------------------|-----|------------------|
| 1 | To instrument rear "pH probe" | 5 | Flow cell holder |
| 2 | Electrode | 6 | Flow cell |
| 3 | Nut | 7 | End cover |
| 4 | Dummy | | |

The flow cell holder can be placed on either the optical unit (for ÄKTApurifier) or on the outside of the valve door (for ÄKTAexplorer). For more information on installation, see ÄKTApurifier User Manual and ÄKTAexplorer Installation Guide.

3.5 Connections

Communication

Connect the network, signal cables and computer according to the electrical drawings in *Electrical and communication connections*, on page 12

Make sure that UNICORN control software is installed on the computer.

Installing Controller unit CU-950

Hang the CU-950 on the left side of the system by inserting the hooks on the front of CU-950 into the channel on the side of the UV-900.

Connect according to diagram in *Electrical and communication connections, on page 12*

When using a fraction collector:

- 1 connect a UniNet-1 cable between Monitor UPC-900/Monitor UV-900 and the fraction collector.
 - 2 connect a termination plug to the empty UniNet-1 socket (Frac-950 only).
-

Flow path

All parts and tubing are mounted on the ÄKTAexplorer and ÄKTAmicro systems at delivery. ÄKTApurifier has no preconnected tubing. It is recommended that the mounting is done by GE Healthcare service engineers. For more information regarding installation, see ÄKTApurifier User Manual.

Setting the delay volume in UNICORN

The delay volume is volume of liquid in the flow path from the UV sensor that identifies the peak to the fraction collector. The length of tubes affects the delay volume that needs to be changed in UNICORN.

To change the delay volume in UNICORN:

- 1 Select **System:Settings** in **System Control**.
- 2 Select **Specials** and then **FracParameters**.
- 3 Enter the delay volume and click **OK**.

Note: *To prevent bacterial growth, the system flow path is filled with 20% ethanol at delivery.*

- Connect tubes for reagents, solvents and sample collection to the correct inlet and outlet connections on the system. For more information, see ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals.
-

Electrical power

Connect the power cord to a grounded power outlet specified in *Section 3.1 Site requirements, on page 32*.

3.6 Spare parts and accessories

For correct up to date information on spare parts and accessories visit:
www.gelifesciences.com/AKTA

4 Operation

About this chapter

This chapter provides instructions for the use of ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro.

4.1 Operation overview

Workflow

The typical workflow in ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro, after turning on the system and connecting it to UNICORN, can be divided into a number of steps.

| Step | Action | Section |
|------|---|---|
| 1 | Create a method | <i>Create a method, on page 43</i> |
| 2 | Prepare the system for a run | <i>Section 4.5 Setting up a run, on page 43</i> <i>Section 4.4 Preparations before start, on page 40</i> |
| 3 | Start a run using a method | <i>Section 4.6 Performing a run, on page 47</i> |
| 4 | During a run - view and change parameters | <i>Viewing the run, on page 48</i> |
| 5 | Procedures after a run | <i>Section 4.7 Procedures after a run, on page 50</i> |
| 6 | Evaluate the results | See UNICORN user documentation. |

Liquid flow path

See *Appendix A Connection diagram - Liquid flow path, on page 74* for an illustration of the liquid flow path in ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro.

4.2 Starting the instrument

Make sure that the waste container and needed buffer bottles are correctly connected. Check that all tubing connections are properly tightened and that all valves are connected to a tube or termination.

Turn on the **Power** switch on the instrument.

4.3 Starting the control system

Starting UNICORN

- 1 Turn on the monitor, computer and optional printer according to the manufacturer's instructions. Wait for the computer to start up.
- 2 Verify that the power indicator on the CU-950 is on when the computer has been turned on.
- 3 Log on to Windows operating system.
- 4 Start UNICORN by double-clicking on the UNICORN shortcut icon on the desktop.



- 5 In the **Logon** dialog, select a user from the **User name** list and enter the password. If you log on for the very first time, select user **default** and enter the password **default**. Click **OK**.




UNICORN starts and the **UNICORN Manager** window opens, see *Figure 4.1*.

Note: See the *UNICORN user documentation* for instructions about how to create new users.

4 Operation

4.3 Starting the control system



| No. | Description |
|-----|--|
| 1 | The Instant Run icon immediately starts the system control wizard used to start a run. |
| 2 | The New Method icon opens the Method Editor module and displays the New Method dialog box. |
| 3 | The System Control icon activates the System Control module and displays the Manual instruction dialog box. |

Figure 4.1: The **UNICORN Manager** window.

Control system in UNICORN

To open the **System Control** module in UNICORN, click the **System Control** icon in the **UNICORN Manager** window, see Figure 4.1.

4.4 Preparations before start

Prepare buffers, solutions and inlets

- 1 Prepare buffers and solutions required for the run.


- 2 Immerse all inlet tubing in the appropriate liquid containers as described in the method.
-

Purging the pump and inlet tubing


Fill the pump and inlet tubing with liquid if small amounts of air need to be removed or if the inlet tubing is empty.

To fill the inlet tubing manually in **System Control**:

- 1 Make sure that no method has been started.
- 2 Set a low flow in **System Control:Manual:Pump:Flow**, for example 0.5 ml/min.
- 3 Click **Execute**.
- 4 Set the inlet valve to the appropriate position in **System Control:Manual:Flowpath:InletValve**, for example inlet **A11** to **A18** or **A2**.
- 5 Connect a syringe to the purge valve.




- 6 Turn the purge valve counter clockwise half a turn to open it.



4 Operation

4.4 Preparations before start

- 7 Slowly draw solution into the syringe. When fluid starts to enter the syringe, continue to draw a few milliliters before closing the purge valve. Check that there is no visible air left in the tubing.



- 8 Repeat for the other purge valve.
 - 9 To fill inlet **B1** and **B2**:
 - a In **System Control:Manual:Pump:Gradient**, select **Target 100% B** and inlet **B1** to fill **B1** or inlet **B2** to fill **B2**. Wait for the valve to turn (a clicking sound) before starting the purging procedure.
 - b When all inlets are filled, click **End**.
-

Connect columns and Superloop™

For column positions, see the method.

Remove air before connecting columns


Air remaining in the system may be removed by purging the pump and by selecting **Pump Wash** and **System Wash**.

- 1 Immerse **A1** tubing in the buffer to be used.
- 2 Select **System Control:Manual:Pump:PumpWash**.

Connecting tubing to columns

Refer to column manufacturer's instructions.

Column attachment drop-to-drop



Attach columns manually by starting a low flow (see *Purging the pump and inlet tubing, on page 41*) and selecting **System Control:Manual:Flowpath:ColumnPosition**.

Preparing the fraction collector

Place the rack chosen in the method in the fraction collector and fill it with appropriate tubes and/or deep well plates.


4.5 Setting up a run

Create a method

To create a method:


- 1 Click the **New Method** icon in the **UNICORN Manager** window, see *Figure 4.1*.

The **New Method** dialog opens.



- 2 If several systems are available, select which system to use in the **For system** list box.
- 3 Select **Wizard** to create a method using the **Method Wizard**. Click **OK**.

The **Method Wizard** dialog opens.



- 4 In the **Method Wizard** dialog:
 - Select the appropriate chromatographic technique in the **Main Selection** list box, for example **Affinity**.
 - Select the column you intend to use in the **Column** list box. The correct column volume, the recommended flow rate, and the correct pressure limit for that column will then be automatically implemented in the method.
 - Select column position in the **Column Position** list box.
 - If required, select **Flexible Flow Rates** and/or **Flow Regulation of the System Pump** and/or **BufferPrep**.
- 5 Click **Next**.
- 6 On each new page in the **Method Wizard**, select the appropriate parameters and click **Next** to continue.

- 7 On the last page, click **Finish**.

The **Run Setup** dialog opens with the **Variables** tab selected by default.

| Block | | Variable | Value | Range |
|---------------------------|----------------------------|-----------|-----------------------------------|-------------------|
| Main | | Column | HiTrap_Protein_A_HP_5_ml [Global] | |
| Flow_Rate | Flow_Rate (ml/min) | 5.000 | | 0.000 - 10.000 |
| Column_Pressure_Limit | Column_PressureLimit (MPa) | 0.30 | | 0.00 - 25.00 |
| Start_Instructions | Wavelength_1 (nm) | 280 | | 190 - 700 |
| | Wavelength_2 (nm) | OFF | | 190 - 700 |
| | Wavelength_3 (nm) | OFF | | 190 - 700 |
| Alarm_Sample_PressLimit | Sample_PressLimit (MPa) | 0.30 | | 0.00 - 1.00 |
| BufferValve_A1_Inlet | BufferValve_A1_Inlet | A11 | | |
| Eluent_A_Inlet | Pump_A_Inlet | A1 | | |
| Eluent_B_Inlet | Pump_B_Inlet | B1 | | |
| Column_Valve | Column_Position | Position2 | | |
| Flush_Sample_ | Sample_FlushVolume (ml) | 0.00 | | 0.00 - 999999.00 |
| Sample_Inlet | Sample_Inlet | S1 | | |
| Flowthrough_Fractionation | Flowthrough_FracSize (ml) | 0.000 | | 0.000 - 99999.000 |
| Direct_SampleLoading | VolumeOfSample (ml) | 5.0 | | 0.0 - 20000.0 |
| Fractionation_Segment_1 | 1_Fraction_Size (ml) | 0.000 | | 0.000 - 99999.000 |
| | 1_PeakFraction_Size (ml) | 0.000 | | 0.000 - 99999.000 |

Show details
 Show unused variables
 Display tooltip for extended variable cells

Edit Variable... Help

- 8 The method is represented by a number of blocks on the **Variables** tab. The blocks are typical steps in a chromatographic run.

Each block contains a number of method variables. If necessary, change the variables to suit your application.

- 9 In the **Method Editor**, select **File:Save As** to save the method.

The **Save As** dialog opens.

Save As

c:\...\Default

| Name | System | Size | Type |
|------|--------|------|-------------|
| | | | Prev Folder |

Method name: Untitled

For System:

- AA_E10S
- AA_PILOT
- AL_E100
- KB-E10
- KB-E10S
- KS_E100
- MH_E100
- RL_E10xT2

Technique:

- Any
- Anion_Exchange
- Cation_Exchange
- RPC
- HIC
- Size_Exclusion
- Affinity
- Chromatofocusing

OK Cancel Help

- 10 In the **Save As** dialog:
 - Enter method name in the **Method name** field and select folder to save the method in.
 - If you have more than one system connected to the computer, select system in the **For System** area for which the method is intended.
 - Select technique in the **Technique** area for which the method was written.
 - 11 Click **OK**.

The method is saved. It can now be started from the **System Control** module.
-

4.6 Performing a run

- 1 **Select method**
 - a In **System Control**, select **File:Run**.
 - b Select the required **Method** from the list.
 - 2 **Specify variables**

Enter identification names for the samples via the keyboard.
 - 3 **Edit result file location and names**

If required, edit the folder path and file names of the result files to be created.
 - 4 **Preparations completed?**

Make sure that the preparations according to *Section 4.4 Preparations before start, on page 40* has been performed.
 - 5 **Check the flow path**

Make sure that:

 - there is enough buffer available
 - the correct inlet is placed in each buffer
 - the outlets are placed in correct bottles
 - the columns are placed in correct positions
 - the chosen fraction collector rack is filled and is in correct place.
 - 6 **Prepare the samples**

The samples should have been prepared and clarified using centrifugation and/or filtration through a 0.45 µm filter¹.

 - a Place the sample in chosen liquid container, vials for the autosampler or fill the capillary loop or the Superloop with sample depending on chosen method.
- 1 If using HisTrap™ FF crude, clarification is not needed.

4 Operation

4.6 Performing a run

Make sure that no air enters the tubing. Place the tubing close to the bottom of the liquid container but not too tight against the bottom.

b Secure the tubing.

7 Final check


Perform a final check that tubing, columns, solutions and buffers are placed according to the method.

8 Start the run

Click **START** to start the run on the selected systems.

Viewing the run

The progress of the run can be viewed in detail in the **System Control** module.



Up to four view panes, **Run Data**, **Curves**, **Flow scheme** and **Logbook** can be displayed showing different aspects of the run in real-time.

- The **Run Data** view pane displays the current values for selected run parameters.

- The **Curves** view pane displays the monitor signal values graphically.
- The **Flow scheme** view pane displays a graphical representation of the chromatography system that shows the current status of the run. During a run, the flow scheme shows open flow path(s) in color and monitor signals with numerical displays.
- The **Logbook** view pane shows the actions as the run proceeds. All actions and unexpected conditions are logged, with date, time and current user name. The log book provides a complete history of the run and is saved in the result file.

Customize the view panes

To customize the view panes, right-click in the respective view pane and select **Properties**. For more information about customizing the view panes, see the UNICORN user documentation.

Ending the run

To stop the run on a system before it is finished:

Click **End** above the **Run data** view pane.

Status indicator colors

The status indicator is located at the bottom of **System Control**.

The table below shows how the indicator colors relate to the run status.

| Indicator color | Run status |
|-----------------|---------------|
| White | End |
| Green | Run or Manual |
| Yellow | Hold |
| Red | Pause |

Error indication

When a warning or an alarm is issued from a system, an error code is displayed. See ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals for guidance.

Evaluate the results

See ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals and UNICORN user documentation for how to evaluate the results.

4.7 Procedures after a run

Cleaning system

To keep the system in good shape, it is important to clean both the tubing and the outside of the system regularly.

- 1 In the **Method Editor** module in UNICORN, create a method for cleaning the system.
 - 2 Wash the outside of the inlet tubings with water and/or ethanol.
 - 3 Immerse the tubing ends to be used in the container with cleaning solution.
 - 4 If the column valve is to be cleaned, remove the columns and reconnect the tubings to the column valves.
 - 5 Run the cleaning method as described in *Section 4.6 Performing a run, on page 47*.
-

Cleaning columns

When running different types of purification methods and different samples after each other, the columns should be cleaned between the runs according to the column instructions. This will remove unspecific bound proteins and prevent column clogging.

- 1 In the **Method Editor** module in UNICORN, create a method for column cleaning in place (CIP).
 - 2 Immerse the tubing ends to be used in the correct containers according to the method for the chosen run.
 - 3 Run the cleaning in place method as described in *Section 4.6 Performing a run, on page 47*.
-

5 Maintenance

About this chapter

This chapter provides instructions for routine component maintenance and a maintenance schedule.

5.1 General

Regular maintenance is important for safe and trouble-free operation of your instrument. The user should perform daily and monthly maintenance. Preventive maintenance should be performed on a yearly basis by qualified service personnel.

For maintenance of a specific component, carefully read the component manual and follow the instructions.



WARNING

Electrical shock hazard. All repairs should be done by service personnel authorized by GE Healthcare. Do not open any covers or replace parts unless specifically stated in the user documentation.



WARNING

Disconnect power. Always disconnect power from the instrument before replacing any component on the instrument, unless stated otherwise in the user documentation.



WARNING

Hazardous chemicals during maintenance. When using hazardous chemicals for system or column cleaning, wash the system or columns with a neutral solution in the last phase or step.

**WARNING**

Do not perform any type of maintenance work while the system is powered electrically or when the piping system is pressurized. Note that the piping system can be pressurized even when the system is closed down.

**WARNING**

When using hazardous chemical and biological agents, take all suitable protective measures, such as wearing protective glasses and gloves resistant to the substances used. Follow local and/or national regulations for safe operation and maintenance of ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro.

**NOTICE**

Cleaning. Keep the instrument dry and clean. Wipe regularly with a soft damp tissue and, if necessary, a mild cleaning agent. Let the instrument dry completely before use.

5.2 User maintenance schedule

Table 5.1 provides a guide to maintenance operations and intervals at which these operations should be performed by the user. The user is however responsible for deciding the type of operations and length of intervals necessary to maintain system function and safety.

Table 5.1: User maintenance schedule

| Interval | Action | Instructions/reference |
|----------|-----------------------------------|---|
| Daily | Leak inspection | Visually inspect the system for leaks. |
| | Wash the system flow path | <ol style="list-style-type: none"> 1 For cleaning the flow path, see <i>Cleaning-In-Place</i>, on page 55. 2 For leaving the system for a few days, see <i>Section 5.8 Storage</i>, on page 58. |
| | Calibrate pH electrode (optional) | Calibrate the pH electrode (if applicable) according to Monitor pH/C-900 User Manual. |

| Interval | Action | Instructions/reference |
|----------|--|---|
| Weekly | Check inlet filters | Check the inlet filters visually and replace them if necessary. |
| | Replace on-line filter (if applicable) | Replace the on-line filter. |
| | Change pump rinsing solution | <p>Change rinsing solution. Always use 20% ethanol with 10 mM NaOH as rinsing solution.</p> <p>If the volume of rinsing solution in the storage bottle has increased, it can be an indication of internal pump leakage. Replace the piston seals according to the User manual.</p> <p>If the volume of rinsing solution in the storage bottle has decreased significantly, check if the rinsing system connectors are mounted properly.</p> <p>If the rinsing system connectors are not leaking, the rinsing membranes or piston seals may be leaking. Replace the membranes and piston seals according to the User manual.</p> |

5 Maintenance

5.2 User maintenance schedule

| Interval | Action | Instructions/reference |
|----------|------------------|---|
| Monthly | Flow restrictor | <p>Check that flow restrictor generates the following back-pressure:</p> <p>FR-904: 0.4 ±0.05 MPa</p> <p>Check the back-pressure as follows:</p> <ol style="list-style-type: none">1 Disconnect the flow restrictor.2 Connect a tubing (approx. 1 m, i.d. 1 mm) to a free port in the injection valve. Set the valve manually to this port. Put the open end in a waste container.3 Run the pump at 10 ml/min with water. Note the back-pressure (Bp1) on the pump display, or in the Run Data window.4 Connect the flow restrictor to the open end of the tubing (observe the IN marking). Put the flow restrictor in the waste container.5 Run the pump at 10 ml/min with water. Note the back-pressure (Bp2) on the pump display, or in the Run Data window.6 Calculate the back-pressure generated by the flow restrictor. Replace it if it is not within limit. |
| Yearly | Valve inspection | <p>Check for external or internal leakage. Replace channel plate and distribution plate yearly or when required. Refer to the relevant valve instruction sheet.</p> |

5.3 Cleaning

Cleaning before planned maintenance/service

To ensure the protection and safety of service personnel, all equipment and work areas must be clean and free of any hazardous contaminants before a Service Engineer starts maintenance work.

Please complete the checklist in the *On Site Service Health & Safety Declaration Form* or the *Health & Safety Declaration Form for Product Return or Servicing*, depending on whether the instrument is going to be serviced on site or returned for service, respectively.

Copy the form you need from *Section 7.4 Health and Safety Declaration Form, on page 71* or print it from the PDF file available on the User Documentation CD.

Cleaning-In-Place

All components in the system are designed for ease of CIP.

After repeated separation cycles, contaminating material might progressively build up in the system and on the column. This material may not have been removed by the cleaning step described above. The nature and degree of contamination depends on the sample and the chromatographic conditions employed. These should be considered when designing a cleaning protocol.

Routine cleaning should be performed at intervals aimed at prevention rather than cleaning the system from growth or contamination.



WARNING

Make sure that the piping system is completely leakage free before performing any CIP on the system.

Make sure that the process control method for cleaning flushes all possible flow paths in the system. After cleaning, rinse the entire system with water or suitable liquid until the piping/tubing system is completely free from the CIP solution (monitors in the system can be used as detectors). Do not leave NaOH or other cleaning agents in the system for long periods.



WARNING

Hazardous chemicals during maintenance. When using hazardous chemicals for system or column cleaning, wash the system or columns with a neutral solution in the last phase or step.



WARNING

NaOH is corrosive and therefore dangerous to health. When using hazardous chemicals, avoid spillage and wear protective glasses and other suitable Personal Protective Equipment (PPE).

See also *Section 5.8 Storage, on page 58*.

5.4 Component maintenance

Maintenance and preventive replacement of parts of the major components are described in the respective manuals included in the system documentation.

The system documentation also includes a spare part list to be used to find common spare parts and their code numbers for ordering. This list can also be found online at www.gelifesciences.com/AKTA.

5.5 Disassembly and assembly of components and consumables

The operator must carefully read and understand the instructions supplied for each component before disassembly and assembly of the component. When replacing consumables, such as tubing and tubing connectors, all necessary safety precautions must be taken. Contact your local GE Healthcare representative if further information or help is needed.



WARNING

Disconnect power. Always disconnect power from the instrument before replacing any component on the instrument, unless stated otherwise in the user documentation.



WARNING

Before disassembly, check that there is no pressure in the piping system.



WARNING

After assembly, the piping system must be tested for leakage at maximum pressure for continued protection against injury risks due to fluid jets, burst pipes or explosive atmosphere.

5.6 Replacement of fuses



WARNING

Disconnect power. Always disconnect power from the instrument before replacing fuses.

Refer to *Section 7.1 Specifications, on page 66* for information about the fuse type and rating. If a fuse repeatedly blows, switch off the system mains switch and contact your local GE Healthcare representative.



WARNING

For continued protection from fire hazard, replace only with same type and rating of fuse.

5.7 Calibration

The table below lists the type and frequency of calibrations that can be done on the instrument. Refer to UNICORN user documentation and to the individual component User Manuals and Instructions for descriptions of how to perform these calibrations. The calibrations are performed from UNICORN by selecting **System:Calibrate** in **System Control**.

| Component | | How often |
|----------------------------|------------------------------|---|
| pH monitor (if applicable) | | Every day. |
| Pump (if applicable) | | Whenever the running conditions are changed, e.g. viscosity of sample or buffer, temperature, backpressure etc. If the sample pump is not used frequently it should be calibrated before use. |
| Pressure reading | | When required. |
| Conductivity flow cell | Cell constant | Only necessary if specific conductivity with high accuracy is measured (Cond_Calib). |
| | Temperature | Must be done when changing the conductivity flow cell (Temp). |
| | Entering a new cell constant | Must be done when changing the conductivity flow cell (Cond_Cell). |

5.8 Storage

General recommendation

For storage, the system must first be cleaned as described in *Cleaning-In-Place*, on page 55. After cleaning, the system must be filled with 0.01 M NaOH or 20% ethanol solution.

Columns and media shall be stored according to their respective instructions.

Storage conditions

The following conditions shall be maintained while the system is in storage:

- Temperature: 2°C to 30°C (preferably room temperature)
- Relative humidity: 0% to 95%, non-condensing (preferably low humidity).

After storage, clean the system, calibrate all monitors, and perform a leakage test before using the system.

6 Troubleshooting

6.1 UV curve problems

| Error symptom | Possible cause | Corrective action |
|--|---|--|
| Ghost peak | Dirt or residues in the flow path from previous runs. Air in the eluents. | Clean the system. Make sure air is removed. |
| | Residue in the column from previous runs | Clean the column according to the column instructions. |
| | Incorrect mixer function | Check the mixer function by placing a stirrer bar on top of the mixer housing. The stirrer bar should rotate when the system is in Run mode. The mixer function can also be checked by running the installation test. |
| Noisy UV-signal, signal drift or instability | Bad UV fiber connections | Check the connections of the UV cell optical fiber. Replace if necessary. |
| | Dirty UV cell | Clean the UV cell by flushing Decon™ 90, Deconex™ 11 or equivalent. |
| | Impure buffer | Check if the signal is still noisy with water. |
| | Air in the pump or in the UV cell | Purge the pump according to <i>Pump User Manual</i> . Run a system wash with buffer. |

6 Troubleshooting

6.1 UV curve problems

| Error symptom | Possible cause | Corrective action |
|-----------------|--|--|
| Low sensitivity | Aging UV lamp | Check the lamp run time according to and replace if necessary. Refer to ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals. |
| | UV lamp in wrong position | Check that the lamp position and the filter position are both set to the wavelength to be used, 280 nm or 254 nm. Refer to ÄKTApurifier and ÄKTAmicro User Manuals. Does not apply to ÄKTAexplorer and ÄKTApurifier without UPC. |
| | The theoretical extinction coefficient too low | Calculate the theoretical extinction coefficient of the protein. If it is zero or very low at 280 nm, the protein cannot be detected. |

6.2 Conductivity curve problems

| Error symptom | Possible cause | Corrective action |
|--------------------------------|----------------------------------|--|
| Baseline drift or noisy signal | Air in the pump or the flow cell | Check the flow restrictor after the flow cell. |
| | Leaking tube connections | Tighten the connectors. If necessary, replace the connectors. |
| | Incorrect mixer function | Check the mixer function by placing a stirrer bar on top of the mixer housing. The stirrer bar should rotate when the system is in Run mode. The mixer function can also be checked by running the installation test. |
| | Dirty conductivity cell | Clean the conductivity cell by flushing 1 M NaOH or 20% ethanol. |
| | Column not equilibrated | Equilibrate the column. If necessary, clean the column. Refer to ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals. |

| Error symptom | Possible cause | Corrective action |
|---|------------------------------------|--|
| Conductivity measurement with the same buffer appears to decrease over time | Dirty flow cell | Clean the flow cell according to procedure in <i>Monitor User Manual</i> . |
| | Decrease in ambient temperature | Use a temperature compensation factor. See <i>Monitor User Manual</i> . |
| Waves on the gradient | Incorrect pump function | Check that the pump is operating and is programmed correctly. |
| | Dirty mixing chamber | Check that the mixing chamber is free from dirt or particles. |
| | Insufficient mixing chamber volume | Change to a larger mixing chamber volume if necessary. |
| | Incorrect motor function | Check the motor operation. Place a hand on the mixer and start it by starting the pump at a low flow rate. You should both hear and feel the mixer motor and stirrer when they are spinning. |
| Ghost peaks appear in the gradient profile | Air in the flow cell | Check for loose tubing connections. Use the flow restrictor. |
| Unlinear gradients or slow response to %B changes | Dirty tubing | Wash the tubing and check pump is operating properly. |
| | Incorrect mixer volume | Change to smaller mixer volume. |

6 Troubleshooting

6.2 Conductivity curve problems

| Error symptom | Possible cause | Corrective action |
|-------------------------------|---|---|
| Incorrect or unstable reading | Loose connection of conductivity flow cable | Check that the conductivity flow cell cable is connected properly. |
| | Incorrect pump and valves function | Check that the pump and valves operate correctly. |
| | Incorrect temperature compensation factor | If temperature compensation is being used, check that the temperature sensor is calibrated, and that the correct temperature compensation factor is in use. |
| | Dirty or incorrectly equilibrated column | Check that the column is equilibrated. If necessary clean the column. |
| | Incorrect mixer function | Check the operation of the mixer. The mixer function is checked by placing a stirrer bar on top of the mixer housing. The stirrer bar should rotate when the system is in Run mode. The mixer function can also be checked by running the installation test. |

6.3 pH curve problems

| Error symptom | Possible cause | Corrective action |
|---------------------------|----------------------------------|--|
| No response to pH changes | Faulty electrode connection | Check that the electrode cable is connected properly. |
| | Damaged electrode | The electrode glass membrane may be cracked. Replace the electrode. |
| | Incorrectly connected pH monitor | Check that the pH monitor is correctly connected according to the ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals. |

| Error symptom | Possible cause | Corrective action |
|--|---------------------------------------|---|
| Small response to pH changes | Dirty pH electrode | Clean the pH electrode as detailed in Monitor pH/C-900 User Manual or UPC-900 User Manual. If the problem remains, replace the pH electrode. |
| Slow pH response or Calibration impossible | Contaminated electrode glass membrane | Check the electrode glass membrane. If it is contaminated, clean the electrode following the instructions in Monitor pH/C-900 User Manual or UPC-900 User Manual. |
| | Membrane has dried out | If the membrane has dried out, the electrode may be restored by soaking it in buffer overnight. |

6 Troubleshooting

6.3 pH curve problems

| Error symptom | Possible cause | Corrective action |
|----------------------------------|-----------------------------------|---|
| Incorrect or unstable pH reading | Problem with electrode | <p>Check that the electrode cable is connected properly.</p> <p>Check that the electrode is correctly inserted in the flow cell and, if necessary, hand-tighten the nut.</p> <p>Check that the pH electrode is not broken.</p> <p>Calibrate the pH electrode.</p> <p>Clean the pH electrode if required, see Monitor pH/C-900 User Manual or UPC-900 User Manual.</p> <p>Compare the response of the pH electrode with that of another pH electrode. If the response differ greatly, the electrode may require cleaning or replacement.</p> <p>In organic solvents such as ethanol, methanol and acetonitrile, stable pH measurements are not possible since dehydration of the membrane will occur. It is recommended that the pH electrode is not used in applications using organic solvents. Mount the dummy electrode instead.</p> |
| | Incorrect pump or valve operation | Check that the pump and valves operate correctly. |
| | Air in the flow cell | If air in the flow cell is suspected, tap the flow cell carefully or tilt it to remove the air. Alternatively, flush the cell with buffer at 20 ml/min (ÄKTAexplorer 100 and ÄKTApurifier 100) or 10 ml/min (ÄKTAexplorer 10 and ÄKTApurifier 10) or 0.5 ml/min (ÄKTAmicro) for 1/2 min. Use the flow restrictor FR-902 after the pH electrode. |
| | Static interference | |

| Error symptom | Possible cause | Corrective action |
|--|----------------------------|--|
| | | There may be interference from static fields. Connect the pH flow cell and the rear panel of the monitor using a standard laboratory 4 mm "banana plug" cable. |
| pH values vary with varied back pressure | Problem with the electrode | Replace the pH electrode. |

6.4 Pressure curve problems

| Error symptom | Possible cause | Corrective action |
|---|--|--|
| Erratic flow, noisy baseline signal, irregular pressure trace | Air bubbles passing through or trapped in the pump | Check all connections for leaks. Check that there is sufficient eluent present in the reservoirs. Use degassed solutions. Purge the pump. Follow the instructions in Pump P-900 User Manual. |
| | Inlet or outlet check valves not functioning correctly | Clean the valves according to Pump P-900 User Manual. |
| | Piston seal leaking | Replace the piston seal according to the instructions in Pump P-900 User Manual. |
| | Blockage or part blockage of flow path | Flush through to clear blockage. If necessary, replace tubing. Check inlet tubing filter. It can become clogged if unfiltered buffers or samples are applied. See instructions for flushing through at the end of the run in Pump P-900 User Manual. |

7 Reference information

About this chapter

This chapter contains technical data, regulatory and other information.

7.1 Specifications

| Parameter | ÄKTAexplorer | ÄKTAmicro | ÄKTApurifier |
|--|----------------------------------|----------------------------------|----------------------------------|
| Ingression protection | IP20 | IP20 | IP20 |
| Supply Voltage | 100-120/220-240 V ~, 50 to 60 Hz | 100-120/220-240 V ~, 50 to 60 Hz | 100-120/220-240 V ~, 50 to 60 Hz |
| Power consumption | 600 VA | 370 VA | 600 VA |
| Fuse specification | T 6.3 AL 250 V | T 6.3 AL 250 V | T 6.3 AL 250 V |
| Dimensions (H x W x D) | 450 x 480 x 610 mm | 450 x 480 x 610 mm | 450 x 490 x 610 mm |
| Weight | 66.8 kg | 55 kg | 41 kg |
| Ambient temperature | 4° to 40 °C | 4° to 40°C | 4° to 40°C |
| Relative humidity tolerance (non-condensing) | 10% to 95% | 10% to 95% | 10% to 95% |
| Atmospheric pressure | 84 to 106 kPa (840 to 1060 mbar) | 84 to 106 kPa (840 to 1060 mbar) | 84 to 106 kPa (840 to 1060 mbar) |

7.2 Chemical resistance

| Chemical | Exposure < 1 day | Exposure up to 2 months | CAS no. | EEC no. | Comments |
|--------------|------------------|-------------------------|---------|---------|----------|
| Acetaldehyde | OK | OK | | | |

| Chemical | Exposure < 1 day | Exposure up to 2 months | CAS no. | EEC no. | Comments |
|----------------------|---------------------|-------------------------------|------------|-----------|---|
| Acetic acid, < 5% | OK | OK | | | |
| Acetic acid, 70% | OK | OK | 64-19-7 | 200-580-7 | |
| Acetonitrile | OK | OK | 75-05-8 | 200-835-2 | FFKM, PP and PE swell. |
| Acetone, 10% | OK | Avoid | | | PVDF is affected by long term use. |
| Ammonia, 30% | OK | OK | 7664-41-7 | 231-635-3 | Silicone is affected by long-term use. |
| Ammonium chloride | OK | OK | 12125-02-9 | 235-186-4 | |
| Ammonium bicarbonate | OK | OK | | | |
| Ammonium nitrate | OK | OK | | | |
| Ammonium sulphate | OK | OK | 7783-20-2 | 231-984-1 | |
| 1-Butanol | OK | OK | | | |
| 2-Butanol | OK | OK | | | |
| Citric acid | OK | OK | 29340-81-6 | 249-576-7 | |
| Chloroform | OK | Avoid | | | Kalrez™, CTFE, PP and PE are affected by long term use. |
| Cyclohexane | OK | OK | | | |
| Detergents | OK | OK | | | |
| Dimethyl sulphoxide | Avoid | Avoid | 67-68-5 | 200-664-3 | PVDF is affected by long term use. |
| 1, 4-Dioxane | Avoid | Avoid | | | ETFE, PP, PE and PVDF are affected by long term use. |
| Ethanol, 100% | OK | OK | 75-08-1 | 200-837-3 | |

7 Reference information

7.2 Chemical resistance

| Chemical | Exposure < 1 day | Exposure up to 2 months | CAS no. | EEC no. | Comments |
|----------------------------|---------------------|-------------------------------|-----------|-----------|--|
| Ethyl acetate | OK | Avoid | | | Silicone not resistant. Pressure limit for PEEK decreases. |
| Ethylene glycol, 100% | OK | OK | 107-21-1 | 203-473-3 | |
| Formic acid, 100% | OK | OK | 64-18-6 | 200-579-1 | Silicone not resistant. |
| Glycerol, 100% | OK | OK | 56-81-5 | 200-289-5 | |
| Guanidinium hydrochloride | OK | OK | | | |
| Hexane | OK | Avoid | | | Silicone not resistant. Pressure limit for PEEK decreases. |
| Hydrochloric acid, 0.1 M | OK | OK | 7647-01-0 | 231-595-7 | Silicone not resistant. |
| Hydrochloric acid, > 0.1 M | OK | Avoid | | | Silicone not resistant. Titanium is affected by long term use. |
| Isopropanol, 100% | OK | OK | 67-63-0 | 200-661-7 | |
| Methanol, 100% | OK | OK | 74-93-1 | 200-659-6 | |
| Nitric acid, diluted | OK | Avoid | | | Silicone not resistant. |
| Nitric acid, 30% | Avoid | Avoid | | | Elgiloy™ is affected by long term use. |
| Phosphoric acid, 10% | OK | Avoid | 7664-38-2 | 231-633-2 | Titanium, aluminium oxide and glass are affected by long term use. |
| Potassium carbonate | OK | OK | 584-08-7 | 209-529-3 | |
| Potassium chloride | OK | OK | 7447-40-7 | 231-211-8 | |

| Chemical | Exposure < 1 day | Exposure up to 2 months | CAS no. | EEC no. | Comments |
|--------------------------------------|---------------------|-------------------------------|-----------|-----------|--|
| Pyridine | Avoid | Avoid | | | ETFE, PP and PE not resistant. |
| Sodium acetate | OK | OK | | | |
| Sodium bicarbonate | OK | OK | | | |
| Sodium bisulphate | OK | OK | | | |
| Sodium borate | OK | OK | | | |
| Sodium carbonate | OK | OK | | | |
| Sodium chloride | OK | OK | 7647-14-5 | 231-598-3 | |
| Sodium hydroxide, 2 M | OK | Avoid | 1310-73-2 | 215-185-5 | PVDF and borosilicate glass are affected by long term use. |
| Sodium sulphate | OK | OK | 7757-82-6 | 231-820-9 | |
| Sulphuric acid, diluted | OK | Avoid | | | PEEK and titanium are affected by long term use. |
| Sulphuric acid, medium concentration | Avoid | Avoid | | | |
| Tetrachloroethylene | Avoid | Avoid | | | Silicone, PP and PE are not resistant. |
| Tetrahydrofuran | Avoid | Avoid | | | ETFE, CTFE, PP and PE are not resistant. |
| Toluene | OK | Avoid | | | Pressure limit for PEEK decreases. |
| Trichloroacetic acid, 1% | OK | OK | 76-03-9 | 200-927-2 | |
| Trifluoroacetic acid, 1% | OK | OK | 176-05-1 | 200-929-3 | |
| Urea, 8M | OK | OK | 57-13-6 | 200-315-5 | |

7 Reference information

7.2 Chemical resistance

| Chemical | Exposure < 1 day | Exposure up to 2 months | CAS no. | EEC no. | Comments |
|-----------------------|---------------------|-------------------------------|---------|---------|--|
| o-Xylene and p-Xylene | OK | Avoid | | | PP and PE are affected by long term use. |

7.3 System recommendations

Refer to ÄKTAexplorer, ÄKTApurifier and ÄKTAmicro User Manuals, or contact your local GE Healthcare representative for the most current information.

7.4 Health and Safety Declaration Form

On site service



GE Services

DOC1149542

On Site Service Health & Safety Declaration Form

Service Ticket #:

To ensure the mutual protection and safety of GE Healthcare service personnel and our customers, all equipment and work areas must be clean and free of any hazardous contaminants before a Service Engineer starts a repair. To avoid delays in the servicing of your equipment, please complete this checklist and present it to the Service Engineer upon arrival. Equipment and/or work areas not sufficiently cleaned, accessible and safe for an engineer may lead to delays in servicing the equipment and could be subject to additional charges.

| Yes | No | Please review the actions below and answer "Yes" or "No". Provide explanation for any "No" answers in box below. |
|--|-----------------------|--|
| | | Instrument has been cleaned of hazardous substances. Please rinse tubing or piping, wipe down scanner surfaces, or otherwise ensure removal of any dangerous residue. Ensure the area around the instrument is clean. If radioactivity has been used, please perform a wipe test or other suitable survey. |
| <input type="radio"/> | <input type="radio"/> | Adequate space and clearance is provided to allow safe access for instrument service, repair or installation. In some cases this may require customer to move equipment from normal operating location prior to GE arrival. |
| <input type="radio"/> | <input type="radio"/> | Consumables, such as columns or gels, have been removed or isolated from the instrument and from any area that may impede access to the instrument. |
| <input type="radio"/> | <input type="radio"/> | All buffer / waste vessels are labeled. Excess containers have been removed from the area to provide access. |
| Provide explanation for any "No" answers here: | | |

Equipment type / Product No: Serial No:

I hereby confirm that the equipment specified above has been cleaned to remove any hazardous substances and that the area has been made safe and accessible.

Name in Capital letters:

Company or institution:

Position or job title: Date (Year/month/date): 20... / ... /

Signed:

GE, imagination at work and GE monogram are trademarks of General Electric Company.
GE Healthcare Bio-Sciences Corp, 800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA.
© 2010-12 General Electric Company—All rights reserved. First published April 2010.
28-9800-26 AB 05/2012

Product return



GE Services

DOC1149544

Health & Safety Declaration Form for Product Return or Servicing

Return authorization number: and/or Service Ticket/Request:

To ensure the mutual protection and safety of GE Healthcare personnel, our customers, transportation personnel and our environment, all equipment must be clean and free of any hazardous contaminants before shipping to GE Healthcare. To avoid delays in the processing of your equipment, please complete this checklist and include it with your return.

- 1. Please note that items will NOT be accepted for servicing or return without this form
- 2. Equipment which is not sufficiently cleaned prior to return to GE Healthcare may lead to delays in servicing the equipment and could be subject to additional charges
- 3. Visible contamination will be assumed hazardous and additional cleaning and decontamination charges will be applied

Please specify if the equipment has been in contact with any of the following:

- Yes No Radioactivity (please specify):
- Yes No Infectious or hazardous biological substances (please specify)
- Yes No Other Hazardous Chemicals (please specify)

Equipment must be decontaminated prior to service / return. Please provide a telephone number where GE Healthcare can contact you for additional information concerning the system / equipment.

Telephone No:

Liquid and/or gas in equipment is: Water Ethanol None, empty Argon, Helium, Nitrogen
Liquid Nitrogen Other, please specify:

Equipment type / Product No: Serial No:

I hereby confirm that the equipment specified above has been cleaned to remove any hazardous substances and that the area has been made safe and accessible.

Name in Capital letters:

Company or institution:

Position or job title: Date (Year/month/date): 20 / /

Signed:

To receive a return authorization number or service number, please call local technical support or customer service.

GE, imagination at work and GE monogram are trademarks of General Electric Company. GE Healthcare Bio-Sciences Corp, 800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA. © 2010-12 General Electric Company—All rights reserved. First published April 2010. 28-9800-27 AB 05/2012


7.5 Ordering information

For ordering information visit www.gelifesciences.com/AKTA.

Appendix A Connection diagram - Liquid flow path

Flow path and components


Liquid flow path for ÄKTAexplorer 10



| No. | Description | No. | Description |
|-----|-------------------|-----|----------------------------------|
| 1 | Buffer Valve (V6) | 7 | Column Selection Valve (V2) |
| 2 | Pump A | 8 | Column Selection Valve (V3) |
| 3 | Pump B | 9 | Flow restrictor |
| 4 | Mixer | 10 | Outlet Valve (V4) |
| 5 | On-line filter | 11 | To Fraction Collector (optional) |

| No. | Description | No. | Description |
|-----|----------------------|-----|-------------|
| 6 | Injection Valve (V1) | | |

Liquid flow path for ÄKTAexplorer 100




| No. | Description | No. | Description |
|-----|----------------------|-----|----------------------------------|
| 1 | Buffer Valve (V6) | 9 | Sample Pump |
| 2 | Pump A | 10 | Flow Direction Valve (V7) |
| 3 | Pump B | 11 | Column Selection Valve (V2) |
| 4 | Mixer | 12 | Restrictor |
| 5 | On-line filter | 13 | Column Selection Valve (V3) |
| 6 | Injection Valve (V1) | 14 | Outlet Valve (V4) |
| 7 | Restrictor | 15 | To Fraction Collector (optional) |

A Connection diagram - Liquid flow path


| No. | Description | No. | Description |
|-----|-------------------|-----|-------------|
| 8 | Sample Valve (V5) | | |

Liquid flow path for ÄKTApurifier 10 and 100



| No. | Description | No. | Description |
|-----|----------------------|-----|---|
| 1 | Pump A | 7 | Bypass |
| 2 | Pump B | 8 | Column |
| 3 | Mixer | 9 | Flow restrictor |
| 4 | On-line filter | 10 | To Fraction Collector (optional) |
| 5 | Injection Valve (V1) | 11 | Cond (optional). Connected after the UV cell. |
| 6 | Sample loop | 12 | pH (optional). Connected after the conductivity cell. |

Liquid flow path for ÄKTAmicro



| No. | Description | No. | Description |
|-----|----------------------|-----|----------------------------------|
| 1 | Pump A | 6 | Sample loop |
| 2 | Pump B | 7 | Bypass |
| 3 | Mixer | 8 | Column |
| 4 | On-line filter | 9 | Flow restrictor (supplied) |
| 5 | Injection Valve (V1) | 10 | To Fraction Collector (optional) |

Appendix B Tubing

Tubing specifications for ÄKTAexplorer

Names in the Label column in *Table B.1* and *Table B.2* refer to tubing labels in the liquid flow path connection diagram, see *Liquid flow path for ÄKTAexplorer 10*, on page 74 and *Liquid flow path for ÄKTApurifier 10 and 100*, on page 76, respectively.

Table B.1: Tubing specifications for ÄKTAexplorer 10

| Use | Label | Material | Length (mm) | I.D. (mm) | Volume (µl) |
|------------------------------|--------|----------|-------------|-----------|-------------------|
| Inlets A11 to A18 | A11-18 | FEP | 1250 | 1.6 | 2.5×10^3 |
| Inlet A1 | A1 | FEP | 750 | 1.6 | 1.5×10^3 |
| Inlet A2 | A2 | FEP | 2000 | 1.6 | 4.0×10^3 |
| Inlet B1 | B1 | FEP | 1800 | 1.6 | 3.6×10^3 |
| Inlet B2 | B2 | FEP | 1800 | 1.6 | 3.6×10^3 |
| Switch valve - pump | A3, B3 | FEP | 150 | 1.6 | 302 |
| Pump - Mixer | G1, G2 | PEEK | 300 | 0.50 | 59 |
| Mixer - On-line filter | G3 | PEEK | 150 | 0.50 | 29 |
| Filter - V1 (Inj. Valve) | G4 | PEEK | 460 | 0.50 | 90 |
| V1 - V2 | G5 | PEEK | 270 | 0.25 | 13 |
| V2 - V3 (bypass) | G6 | PEEK | 620 | 0.50 | 122 |
| V3 - UV | G7 | PEEK | 550 | 0.25 | 27 |
| UV - Cond. Cell | G8 | PEEK | 160 | 0.25 | 8 |
| Cond. Cell - Flow restrictor | G9 | PEEK | 450 | 0.25 | 22 |
| Flow restrictor - V4 | G10 | PEEK | 120 | 0.25 | 6 |
| V4 - Fraction collector | G11 | PEEK | 600 | 0.25 | 29 |
| V4 - Flow through | F3 | PEEK | 1000 | 0.50 | 196 |
| V1 - Waste | W1, W2 | ETFE | 1300 | 0.75 | 574 |

| Use | Label | Material | Length (mm) | I.D. (mm) | Volume (µl) |
|------------|-------|----------|-------------|-----------|-------------|
| V4 - Waste | W3 | ETFE | 1000 | 0.75 | 442 |

Table B.2: Tubing specifications for ÄKTAexplorer 100

| Use | Label | Material | Length (mm) | I.D. (mm) | Volume (µl) |
|--|--------|----------|-------------|-----------|--------------------|
| Inlets A11 to A18 | A11-18 | FEP | 1250 | 2.9 | 8.3×10^3 |
| Inlet A1 | A1 | FEP | 750 | 2.9 | 4.9×10^3 |
| Inlet A2 | A2 | FEP | 2000 | 2.9 | 13.2×10^3 |
| Inlet B1 | B1 | FEP | 1800 | 2.9 | 11.9×10^3 |
| Inlet B2 | B2 | FEP | 1800 | 2.9 | 11.9×10^3 |
| Switch valve - pump | A3, B3 | FEP | 150 | 2.9 | 991 |
| Pump - Mixer | G1, G2 | PEEK | 300 | 0.75 | 133 |
| Mixer - On-line filter | G3 | PEEK | 150 | 0.75 | 66 |
| Filter - V1 (Inj. Valve) | G4 | PEEK | 460 | 0.75 | 203 |
| V1 - V7 | G5 | PEEK | 470 | 0.75 | 208 |
| V7 - V2 | G6 | PEEK | 410 | 0.75 | 181 |
| Bypass | G7 | PEEK | 620 | 0.75 | 274 |
| V3 - V7 | G8 | PEEK | 470 | 0.75 | 208 |
| V7 loop | G9 | PEEK | 180 | 0.75 | 80 |
| V7 - UV | G10 | PEEK | 370 | 0.75 | 163 |
| UV - Cond. Cell | G11 | ETFE | 160 | 0.75 | 71 |
| Cond. Cell - pH Cell | G12 | ETFE | 450 | 0.75 | 199 |
| pH Cell - Flow restrictor ¹ | G13 | ETFE | 110 | 0.75 | 49 |
| Flow restrictor - V4 ¹ | G14 | ETFE | 120 | 0.75 | 53 |
| V4 - Fraction collector | G15 | ETFE | 600 | 0.75 | 265 |
| Flow through | F3 | ETFE | 1000 | 1.0 | 785 |
| V1 - Waste | W1 | ETFE | 1300 | 0.75 | 574 |

B Tubing

| Use | Label | Material | Length (mm) | I.D. (mm) | Volume (μl) |
|---------------------|-------|----------|-------------|-----------|-------------------|
| Sample pump - Waste | W2 | FEP | 1300 | 1.6 | 2.6×10^3 |
| V4 - Waste | W3 | ETFE | 1000 | 1.0 | 785 |

¹ Not mounted at factory

Tubing specifications for ÄKTApurifier

Names in the Label column in *Table B.3* and *Table B.4* refer to tubing labels in the liquid flow path connection diagram, see *Liquid flow path for ÄKTApurifier 10 and 100, on page 76*.

Table B.3: Tubing specifications for ÄKTApurifier 10

| Use | Label | Material | Length (mm) | I.D. (mm) | Volume (μl) |
|------------------------------|------------|----------|-------------|-----------|-------------------|
| Inlet A1 | A1 | FEP | 1700 | 1.6 | 3.4×10^3 |
| Inlet A2 | A2 | FEP | 1900 | 1.6 | 3.8×10^3 |
| Inlet B1 | B1 | FEP | 1500 | 1.6 | 3.0×10^3 |
| Inlet B2 | B2 | FEP | 1700 | 1.6 | 3.4×10^3 |
| Switch valve - pump | A3, B3 | FEP | 150 | 1.6 | 302 |
| Pump - Mixer | G1, G2 | PEEK | 330 | 0.50 | 65 |
| Mixer - On-line filter | G3 | PEEK | 200 | 0.50 | 39 |
| Filter - V1 (Inj. Valve) | G4 | PEEK | 180 | 0.50 | 35 |
| V1 - UV (bypass) | G5 | PEEK | 400 | 0.50 | 79 |
| UV - Cond. Cell | G6 | PEEK | 160 | 0.25 | 8 |
| Cond. Cell - Flow restrictor | G7 | PEEK | 80 | 0.25 | 4 |
| Flow restrictor - V4 | G8 | PEEK | 140 | 0.25 | 7 |
| V4 - Fraction collector | G9 | PEEK | 600 | 0.25 | 29 |
| V4 - Flow through | F3 | PEEK | 1000 | 0.50 | 196 |
| V1/V4 - Waste | W1, W2, W3 | ETFE | 100 | 0.75 | 442 |

Table B.4: Tubing specifications for ÄKTApurifier 100

| Use | Label | Material | Length (mm) | I.D. (mm) | Volume (µl) |
|------------------------------|------------|----------|-------------|-----------|--------------------|
| Inlet A1 | A1 | FEP | 1700 | 2.9 | 11.2×10^3 |
| Inlet A2 | A2 | FEP | 1900 | 2.9 | 12.5×10^3 |
| Inlet B1 | B1 | FEP | 1500 | 2.9 | 9.9×10^3 |
| Inlet B2 | B2 | FEP | 1700 | 2.9 | 11.2×10^3 |
| Switch valve - pump | A3, B3 | FEP | 150 | 2.9 | 991 |
| Pump - Mixer | G1, G2 | PEEK | 330 | 0.75 | 146 |
| Mixer - On-line filter | G3 | PEEK | 200 | 0.75 | 88 |
| Filter - V1 (Inj. Valve) | G4 | PEEK | 180 | 0.75 | 80 |
| V1 - UV (bypass) | G5 | PEEK | 400 | 0.75 | 177 |
| UV - Cond. Cell | G6 | ETFE | 160 | 0.75 | 71 |
| Cond. Cell - Flow restrictor | G7 | ETFE | 80 | 0.75 | 35 |
| Flow restrictor - V4 | G8 | ETFE | 140 | 0.75 | 62 |
| V4 - Fraction collector | G9 | ETFE | 600 | 0.75 | 265 |
| V4 - Flow through | F3 | ETFE | 1000 | 1.0 | 785 |
| V1/V4 - Waste | W1, W2, W3 | ETFE | 100 | 0.75 | 442 |

Tubing specifications for ÄKTAmicro

Names in the Tubing i.d. column in *Table* refer to tubing labels in the liquid flow path connection diagram, see *Liquid flow path for ÄKTAmicro*, on page 77.

| Tubing i.d. | Tubing o.d. | Material | Color | Max. pressure | Volume of 10 cm | Connection points |
|---------------------------------|-------------|--------------|--------|---------------|-----------------|--|
| 0.35 mm (G1, G2) | 1.6 mm | Titanium | Grey | > 35 MPa | 9.6 µl | From Pump P-905 to Mixer |
| Union, m/m | – | PEEK | Black | > 35 MPa | – | Between Mixer and On-line filter Between UV cell and Conductivity cell |
| 0.15 mm (G3 to G5) ¹ | 1/16" | PEEK | Violet | > 35 MPa | 1.8 µl | From On-line filter to outlet (or to fraction collector/MS, optional). (PEEK tubing i.d. 0.15 mm is installed at delivery) |
| 0.15 mm (Di. 0.15) | 375 µm | Fused silica | Brown | > 35 MPa | 1.8 µl | From On-line filter to outlet (or to fraction collector/MS, optional). (Tubing kit 0.15 is supplied with the system) |
| 0.10 mm (Di. 0.10) | 200 µm | Fused silica | Brown | > 35 MPa | 0.8 µl | From On-line filter to outlet (or to fraction collector/MS, optional). (Tubing kit 0.10 is supplied with the system) |
| 0.75 mm (W1, W2) | 1/16" | ETFE | Clear | 7 MPa | 44.2 µl | Waste tubing |
| 1.6 mm (A1 to A3, B1 to B3) | 1/8" | FEP | Clear | 3.4 MPa | 201.1 µl | Inlet tubing to Pump P-905 |

¹ 0.25 and 0.5 mm i.d. tubing is supplied for non-analytical applications.

For local office contact information, visit
www.gelifesciences.com/contact

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden

www.gelifesciences.com/AKTA

GE, imagination at work and GE monogram are trademarks of General Electric Company.

ÅKTAexplorer, ÅKTAmicro, ÅKTApurifier, HisTrap, Superloop and UNICORN are trademarks of GE Healthcare companies.

Decon is a trademark of Decon Laboratories Ltd.

Deconex is a trademark of Borer Chemie AG.

Elgiloy is a trademark of Elgiloy Limited Partnership.

Kalrez is a trademark of DuPont Performance Elastomers L.L.C.

Microsoft and Windows are trademarks of Microsoft Corporation.

Any use of UNICORN is subject to GE Healthcare Standard Software End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.

UNICORN © 2003-2013 General Electric Company.

© 2009-2013 General Electric Company – All rights reserved.

First published Jul. 2009

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Limited
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Japan Corporation
Sanken Bldg. 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan



imagination at work